

INCLUSIVE[®]
DENTAL IMPLANT SYSTEM

DETAIL C
SCALE 25:1



Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2003

This is to certify that:

Prismatik Dentalcraft, Inc.
a wholly owned subsidiary of
Glidewell Laboratories
2212 Dupont Drive
Irvine
California
92612
USA

Holds Certificate No:

FM 573328

and operates a Quality Management System which complies with the requirements of ISO 13485:2003 for the following scope:

Design and manufacturing of dental restorative products. Design and development, manufacture, and distribution of dental model scanning and milling systems for dental restorations.

For and on behalf of BSI:



Carlos Pitanga, SVP, System Certification and Compliance

Original Registration Date: 09/02/2011

Effective Date: 07/15/2016

Expiry Date: 08/24/2017



**CMDCAS
Recognized
Registrar**



Page: 1 of 2

...making excellence a habit.™

This certificate remains the property of BSI and shall be returned immediately upon request.
An electronic certificate can be authenticated [online](#). Printed copies can be validated at www.bsigroup.com/ClientDirectory
To be read in conjunction with the scope above or the attached appendix.

Introducing the Inclusive® Dental Implant System

The Inclusive® Dental Implant System is the result of an expansive effort to offer clinicians and their patients a simple, convenient, affordable way to restore missing dentition. Choose between conventional-diameter or small-diameter implant assortments, supported by a wide array of restorative options, and discover a system streamlined for ease of use, yet versatile enough to address almost any clinical indication. From cementable or screw-retained single-tooth replacement to full-arch reconstruction with removable or fixed-removable prostheses, the Inclusive Dental Implant System can provide a predictable foundation for clinical success.

About the Manufacturer

Prismatik Dentalcraft was established in 2006 with the mission of making implant dentistry the standard of care for edentulous patients across the economic spectrum. To realize this goal, we carefully assembled a team of experts with decades of combined experience in the design, engineering, and manufacture of dental implants. With a support staff of highly respected researchers, material scientists, clinical specialists, and dental technicians, Prismatik is dedicated to advancing implant therapies by combining proven treatment protocols with progressive materials, technologies, and techniques.

Expert Personnel



Our team of experts have decades of combined experience in the design and manufacture of dental implants.

State-of-the-Art Equipment



Our Swiss-type lathes and multi-axis milling machines are ideal for implants and prosthetics requiring extreme precision.

Made in the U.S.A.



Our ISO-certified facility in Irvine, Calif. operates under FDA Current Good Manufacturing Practices (CGMPs).

5	Tapered Implant System
6	Technical Specifications
7	Tapered Implants
8	Healing Components
9	Impression Components
11	Stock Abutments
12	Screws & Guide Pins
13	Multi-Unit Abutment System
15	LOCATOR® Abutment System
17	Surgical Kit Overview
19	Drills & Taps
20	Drivers & Accessories
21	Mini Implant System
22	Technical Specifications
23	Mini Implant with O-ball Head
24	Prosthetic Components
25	Surgical Kit Overview
27	Instrumentation
28	Implant Packaging
29	Policies and Warranty
30	Parts Index

*Inclusive® is a registered trademark of Prismatic Dentalcraft, Inc.
LOCATOR is a registered trademark of Zest Anchors, LLC.*

Copyright © 2017, Prismatic Dentalcraft, Inc. Prismatic Dentalcraft, Inc. is not responsible for any damages or other liabilities (including attorney fees) resulting, or claimed to result in whole or in part, from actual or alleged problems arising out of the use of this information. The techniques, procedures and theories presented herein are provided in good faith and believed to be correct as of the date hereof. Any dental professional viewing this presentation must make his or her own decisions about the use of the materials and techniques for specific situations.



No representations as to the completeness or accuracy of this information is given, and no representations or warranties, either expressed or implied, of merchantability, fitness for a particular purpose or of any other nature are made here under with respect to the information or the product to which information refers.



INCLUSIVE®
TAPERED IMPLANT SYSTEM







Technical Specifications



Tapered Implants

Inclusive® Tapered Implants are machined from high-strength titanium alloy and incorporate design features for achieving high primary stability, crestal bone preservation, soft tissue maintenance, maximum bone-to-implant contact, natural root shape, and compatibility with popular prosthetic components and instrumentation. Inclusive Tapered Implants are available in five major diameters, each of which is available in five lengths. Three standard prosthetic platform sizes (3.0 mm, 3.5 mm, 4.5 mm) minimizes the number of system components and simplifies the restorative process.

Description		Code
Inclusive® Tapered Implant		
 70-1070-IMP0035 70-1070-IMP0036 70-1070-IMP0032 70-1070-IMP0033 70-1070-IMP0034	3.2 mmD x 8 mmL x 3.0 mmP	70-1070-IMP0035
	3.2 mmD x 10 mmL x 3.0 mmP	70-1070-IMP0036
	3.2 mmD x 11.5 mmL x 3.0 mmP	70-1070-IMP0032
	3.2 mmD x 13 mmL x 3.0 mmP	70-1070-IMP0033
	3.2 mmD x 16 mmL x 3.0 mmP	70-1070-IMP0034
 70-1070-IMP0005 70-1070-IMP0006 70-1070-IMP0007 70-1070-IMP0008 70-1070-IMP0009	3.7 mmD x 8 mmL x 3.5 mmP	70-1070-IMP0005
	3.7 mmD x 10 mmL x 3.5 mmP	70-1070-IMP0006
	3.7 mmD x 11.5 mmL x 3.5 mmP	70-1070-IMP0007
	3.7 mmD x 13 mmL x 3.5 mmP	70-1070-IMP0008
	3.7 mmD x 16 mmL x 3.5 mmP	70-1070-IMP0009
 70-1070-IMP0026 70-1070-IMP0027 70-1070-IMP0028 70-1070-IMP0029 70-1070-IMP0030	4.2 mmD x 8 mmL x 3.5 mmP	70-1070-IMP0026
	4.2 mmD x 10 mmL x 3.5 mmP	70-1070-IMP0027
	4.2 mmD x 11.5 mmL x 3.5 mmP	70-1070-IMP0028
	4.2 mmD x 13 mmL x 3.5 mmP	70-1070-IMP0029
	4.2 mmD x 16 mmL x 3.5 mmP	70-1070-IMP0030
 70-1070-IMP0010 70-1070-IMP0011 70-1070-IMP0012 70-1070-IMP0013 70-1070-IMP0014	4.7 mmD x 8 mmL x 4.5 mmP	70-1070-IMP0010
	4.7 mmD x 10 mmL x 4.5 mmP	70-1070-IMP0011
	4.7 mmD x 11.5 mmL x 4.5 mmP	70-1070-IMP0012
	4.7 mmD x 13 mmL x 4.5 mmP	70-1070-IMP0013
	4.7 mmD x 16 mmL x 4.5 mmP	70-1070-IMP0014
 70-1070-IMP0015 70-1070-IMP0016 70-1070-IMP0017 70-1070-IMP0018 70-1070-IMP0019	5.2 mmD x 8 mmL x 4.5 mmP	70-1070-IMP0015
	5.2 mmD x 10 mmL x 4.5 mmP	70-1070-IMP0016
	5.2 mmD x 11.5 mmL x 4.5 mmP	70-1070-IMP0017
	5.2 mmD x 13 mmL x 4.5 mmP	70-1070-IMP0018
	5.2 mmD x 16 mmL x 4.5 mmP	70-1070-IMP0019

Healing Components

Description	Implant	Code
-------------	---------	------

Inclusive® Tapered Implant Healing Abutment


			3.0 mmP x 3.2 mmD x 3 mmH	Ø3.2	70-1069-IMP0009
			3.0 mmP x 3.2 mmD x 5 mmH	Ø3.2	70-1069-IMP0010
			3.0 mmP x 3.2 mmD x 7 mmH	Ø3.2	70-1069-IMP0011
			3.0 mmP x 3.7 mmD x 3 mmH	Ø3.2	70-1069-IMP0012
			3.0 mmP x 3.7 mmD x 5 mmH	Ø3.2	70-1069-IMP0013
			3.0 mmP x 3.7 mmD x 7 mmH	Ø3.2	70-1069-IMP0014
			3.5 mmP x 3.5 mmD x 3 mmH	Ø3.7/4.2	70-1069-IMP0001
			3.5 mmP x 3.5 mmD x 5 mmH	Ø3.7/4.2	70-1069-IMP0002
			3.5 mmP x 3.5 mmD x 7 mmH	Ø3.7/4.2	70-1069-IMP0015
			3.5 mmP x 4.7 mmD x 3 mmH	Ø3.7/4.2	70-1069-IMP0003
			3.5 mmP x 4.7 mmD x 5 mmH	Ø3.7/4.2	70-1069-IMP0004
			3.5 mmP x 4.7 mmD x 7 mmH	Ø3.7/4.2	70-1069-IMP0016
			4.5 mmP x 4.5 mmD x 3 mmH	Ø4.7/5.2	70-1069-IMP0005
			4.5 mmP x 4.5 mmD x 5 mmH	Ø4.7/5.2	70-1069-IMP0006
			4.5 mmP x 4.5 mmD x 7 mmH	Ø4.7/5.2	70-1069-IMP0017
			4.5 mmP x 5.7 mmD x 3 mmH	Ø4.7/5.2	70-1069-IMP0007
			4.5 mmP x 5.7 mmD x 5 mmH	Ø4.7/5.2	70-1069-IMP0008
			4.5 mmP x 5.7 mmD x 7 mmH	Ø4.7/5.2	70-1069-IMP0018

(Healing Abutments require a .050" hex instrument for placement - see page 20.)


Healing Components

Description	Implant	Code
-------------	---------	------


Inclusive® Tapered Implant Temporary Abutment



3.0 mmP	Ø3.2	70-1073-PRA0048
Non-Engaging, 3.0 mmP	Ø3.2	70-1073-PRA0047



3.5 mmP	Ø3.7/4.2	70-1073-PRA0022
Non-Engaging, 3.5 mmP	Ø3.7/4.2	70-1073-PRA0023




4.5 mmP	Ø4.7/5.2	70-1073-PRA0024
Non-Engaging, 4.5 mmP	Ø4.7/5.2	70-1073-PRA0025


(Temporary Abutments include Titanium Screw and Guide Pin.)

Impression Components


Inclusive® Tapered Implant Transfer Coping



Closed-Tray, 3.0 mmP x 3.2 mmD	Ø3.2	70-1072-PRC0051
Closed-Tray, 3.0 mmP x 3.7 mmD	Ø3.2	70-1072-PRC0053



Closed-Tray, 3.5 mmP x 3.5 mmD	Ø3.7/4.2	70-1072-PRC0007
Closed-Tray, 3.5 mmP x 4.7 mmD	Ø3.7/4.2	70-1072-PRC0009




Closed-Tray, 4.5 mmP x 4.5 mmD	Ø4.7/5.2	70-1072-PRC0011
Closed-Tray, 4.5 mmP x 5.7 mmD	Ø4.7/5.2	70-1072-PRC0013

(Closed-Tray Transfer Copings include Closed-Tray Screw.)

Impression Components

Description	Implant	Code
-------------	---------	------

Inclusive® Tapered Implant Transfer Coping



70-1072-PRC0052



70-1072-PRC0054

Open-Tray, 3.0 mmP x 3.2 mmD

Ø3.2

70-1072-PRC0052

Open-Tray, 3.0 mmP x 3.7 mmD

Ø3.2

70-1072-PRC0054



70-1072-PRC0006



70-1072-PRC0008

Open-Tray, 3.5 mmP x 3.5 mmD


Ø3.7/4.2

70-1072-PRC0006

Open-Tray, 3.5 mmP x 4.7 mmD

Ø3.7/4.2

70-1072-PRC0008



70-1072-PRC0010



70-1072-PRC0012

Open-Tray, 4.5 mmP x 4.5 mmD

Ø4.7/5.2

70-1072-PRC0010

Open-Tray, 4.5 mmP x 5.7 mmD

Ø4.7/5.2

70-1072-PRC0012

(Open-Tray Transfer Copings include Open-Tray Screw and Guide Pin.)

Inclusive® Tapered Implant Clinical Scanning Abutment



70-1073-PRA0051



70-1073-PRA0042



70-1073-PRA0043

Anterior, 3.0 mmP

Ø3.2

70-1073-PRA0051

Anterior, 3.5 mmP

Ø3.7/4.2

70-1073-PRA0042

Anterior, 4.5 mmP

Ø4.7/5.2

70-1073-PRA0043



70-1073-PRA0050



70-1073-PRA0040



70-1073-PRA0041

Posterior, 3.0 mmP

Ø3.2

70-1073-PRA0050

Posterior, 3.5 mmP

Ø3.7/4.2

70-1073-PRA0040

Posterior, 4.5 mmP

Ø4.7/5.2

70-1073-PRA0041

(Scanning Abutments include captured screw.)

Inclusive® Tapered Implant Laboratory Scanning Abutment



70-1073-PRA0049



70-1073-PRA0033



70-1073-PRA0034

3.0 mmP (3-pack)

Ø3.2

70-1073-PRA0049

3.5 mmP (3-pack)

Ø3.7/4.2

70-1073-PRA0033

4.5 mmP (3-pack)

Ø4.7/5.2

70-1073-PRA0034

(Scanning Abutments include captured screw.)

Inclusive® Tapered Implant Analog



70-1072-PRC0055



70-1072-PRC0015



70-1072-PRC0017

3.0 mmP

Ø3.2

70-1072-PRC0055

3.5 mmP

Ø3.7/4.2

70-1072-PRC0015

4.5 mmP






Ø4.7/5.2

70-1072-PRC0017







Stock Abutments

Description Implant Code

Inclusive® Tapered Implant 15° Titanium Esthetic Abutment

		Anterior, 3.0 mmP x 1.5 mmH	Ø3.2	70-1073-PRA0054
70-1073-PRA0054	70-1073-PRA0059	Anterior, 3.0 mmP x 2.3 mmH	Ø3.2	70-1073-PRA0059
		Anterior, 3.5 mmP	Ø3.7/4.2	70-1073-PRA0009
70-1073-PRA0009	70-1073-PRA0031	Anterior, 3.5 mmP x 1.5 mmH	Ø3.7/4.2	70-1073-PRA0031
		Anterior, 4.5 mmP	Ø4.7/5.2	70-1073-PRA0010
70-1073-PRA0010	70-1073-PRA0062	Anterior, 4.5 mmP x 1.5 mmH	Ø4.7/5.2	70-1073-PRA0062
		Posterior, 4.5 mmP	Ø4.7/5.2	70-1073-PRA0013
70-1073-PRA0013				

Inclusive® Tapered Implant Titanium Esthetic Abutment

		Anterior, 3.0 mmP x 1.5 mmH	Ø3.2	70-1073-PRA0055
70-1073-PRA0055	70-1073-PRA0060	Anterior, 3.0 mmP x 2.3 mmH	Ø3.2	70-1073-PRA0060
		Anterior, 3.5 mmP x 1.5 mmH	Ø3.7/4.2	70-1073-PRA0032
70-1073-PRA0032	70-1073-PRA0011	Anterior, 3.5 mmP	Ø3.7/4.2	70-1073-PRA0011
		Anterior, 4.5 mmP	Ø4.7/5.2	70-1073-PRA0012
70-1073-PRA0012	70-1073-PRA0063	Anterior, 4.5 mmP x 1.5 mmH	Ø4.7/5.2	70-1073-PRA0063
		Posterior, 4.5 mmP	Ø4.7/5.2	70-1073-PRA0014
70-1073-PRA0014	70-1073-PRA0015	Posterior, Wide, 4.5 mmP	Ø4.7/5.2	70-1073-PRA0015

(Esthetic Abutments include Titanium Screw.)

Inclusive® Tapered Implant Titanium Abutment

		3.0 mmP x 4.5 mmH	Ø3.2	70-1073-PRA0052
70-1073-PRA0052	70-1073-PRA0053	3.0 mmP x 6 mmH	Ø3.2	70-1073-PRA0053
		3.5 mmP x 4.5 mmH	Ø3.7/4.2	70-1073-PRA0045
70-1073-PRA0045	70-1073-PRA0006	3.5 mmP x 6 mmH	Ø3.7/4.2	70-1073-PRA0006
		4.5 mmP x 4.5 mmH	Ø4.7/5.2	70-1073-PRA0046
70-1073-PRA0046	70-1073-PRA0008	4.5 mmP x 6 mmH	Ø4.7/5.2	70-1073-PRA0008

(Titanium Screw sold separately; please see page 12.)

Stock Abutments

Description	Implant	Code
-------------	---------	------

Inclusive® Tapered Implant UCLA Plastic Abutment

			3.0 mmP	Ø3.2	70-1073-PRA0058
			3.5 mmP	Ø3.7/4.2	70-1073-PRA0018
			4.5 mmP	Ø4.7/5.2	70-1073-PRA0020
70-1073-PRA0058	70-1073-PRA0018	70-1073-PRA0020			

			Non-Engaging, 3.0 mmP	Ø3.2	70-1073-PRA0057
			Non-Engaging, 3.5 mmP	Ø3.7/4.2	70-1073-PRA0019
			Non-Engaging, 4.5 mmP	Ø4.7/5.2	70-1073-PRA0021
70-1073-PRA0057	70-1073-PRA0019	70-1073-PRA0021			

(UCLA Abutments include Titanium Screw and Guide Pin.)


Screws & Guide Pins

	Titanium Screw 3.0 mmP	Ø3.2	70-1072-PRC0058
70-1072-PRC0058 70-1072-PRC0057	Titanium Screw 3.0 mmP (5-pack)	Ø3.2	70-1072-PRC0057
	Titanium Screw	Ø3.7/4.2/4.7/5.2	70-1072-PRC0001
70-1072-PRC0001 70-1072-PRC0003	Titanium Screw (5-pack)	Ø3.7/4.2/4.7/5.2	70-1072-PRC0003
	Guide Pin 3.0 mmP (5-pack)	Ø3.2	70-1072-PRC0056
70-1072-PRC0056	Guide Pin (5-pack)	Ø3.7/4.2/4.7/5.2	70-1072-PRC0002
70-1072-PRC0002			


Multi-Unit Abutment System

Description	Implant	Code
-------------	---------	------

Inclusive® Tapered Implant Multi-Unit Abutment




3.0 mmP x 1.5 mmH	Ø3.2	70-1074-PRS0040
3.0 mmP x 2.5 mmH	Ø3.2	70-1074-PRS0041
3.0 mmP x 3.5 mmH	Ø3.2	70-1074-PRS0042



3.5 mmP x 1 mmH	Ø3.7/4.2	70-1074-PRS0009
3.5 mmP x 2 mmH	Ø3.7/4.2	70-1074-PRS0010
4.5 mmP x 1 mmH	Ø4.7/5.2	70-1074-PRS0011
4.5 mmP x 2 mmH	Ø4.7/5.2	70-1074-PRS0012
4.5 mmP x 3 mmH	Ø4.7/5.2	70-1074-PRS0013
4.5 mmP x 5 mmH	Ø4.7/5.2	70-1074-PRS0014


(Prosthetic Screw sold separately)

Inclusive® Tapered Implant 17° Multi-Unit Abutment



3.5 mmP x 2 mmH	Ø3.7/4.2	70-1074-PRS0001
3.5 mmP x 3 mmH	Ø3.7/4.2	70-1074-PRS0002
4.5 mmP x 2 mmH	Ø4.7/5.2	70-1074-PRS0003
4.5 mmP x 3 mmH	Ø4.7/5.2	70-1074-PRS0004
4.5 mmP x 5 mmH	Ø4.7/5.2	70-1074-PRS0005

Inclusive® Tapered Implant 30° Multi-Unit Abutment



3.5 mmP x 3 mmH	Ø3.7/4.2	70-1074-PRS0006
4.5 mmP x 3 mmH	Ø4.7/5.2	70-1074-PRS0007
4.5 mmP x 5 mmH	Ø4.7/5.2	70-1074-PRS0008

(Angled Multi-Unit Abutments include Angled Abutment Screw. Prosthetic Screw sold separately; please see page 14.)

Multi-Unit Abutment System

	Description	Code
Inclusive® Multi-Unit		
 70-1074-PRS0015	Temporary Healing Cap with Prosthetic Screw .050" Hex	70-1074-PRS0015
 70-1073-PRA0037	Coping 4.5 mmH, Non-Engaging with Prosthetic Screw .050" Hex	70-1073-PRA0037
 70-1074-PRS0030	Gold/Plastic Coping	70-1074-PRS0030
 70-1074-PRS0018	Gold Bar Coping	70-1074-PRS0018
 70-1074-PRS0052	Impression Coping, Open-Tray with Guide Pin .050" Hex	70-1074-PRS0052
 70-1074-PRS0021	Impression Coping, Closed-Tray	70-1074-PRS0021
 70-1074-PRS0025	Titanium Temporary with Prosthetic Screw .050" Hex	70-1074-PRS0025
 70-1074-PRS0067	Plastic Coping with Prosthetic Screw .050" Hex	70-1074-PRS0067
 70-1074-PRS0016	Abutment Analog	70-1074-PRS0016
 70-1074-PRS0028	Driver	70-1074-PRS0028
 70-1074-PRS0061	Guide Pin .050" Hex	70-1074-PRS0061
 70-1072-PRC0014	Prosthetic Screw .050" Hex	70-1072-PRC0014
 70-1074-PRS0031	Inclusive® Angled Multi-Unit Abutment Screw	70-1074-PRS0031

LOCATOR® Abutment System

	Description	Implant	Code
LOCATOR® Abutment			
 70-1075-ZST0065	3.0 mmP x 1 mmH Cuff [#2212]	Ø3.2	70-1075-ZST0065
 70-1075-ZST0066	3.0 mmP x 3 mmH Cuff [#2214]	Ø3.2	70-1075-ZST0066
 70-1075-ZST0011	3.5 mmP x 0 mmH Cuff [#8661]	Ø3.7/4.2	70-1075-ZST0011
 70-1075-ZST0012	3.5 mmP x 1 mmH Cuff [#8662]	Ø3.7/4.2	70-1075-ZST0012
 70-1075-ZST0013	3.5 mmP x 2.5 mmH Cuff [#8663]	Ø3.7/4.2	70-1075-ZST0013
 70-1075-ZST0014	3.5 mmP x 3.5 mmH Cuff [#8664]	Ø3.7/4.2	70-1075-ZST0014
 70-1075-ZST0015	4.5 mmP x 0 mmH Cuff [#8671]	Ø4.7/5.2	70-1075-ZST0015
 70-1075-ZST0016	4.5 mmP x 1 mmH Cuff [#8672]	Ø4.7/5.2	70-1075-ZST0016
 70-1075-ZST0017	4.5 mmP x 2 mmH Cuff [#8673]	Ø4.7/5.2	70-1075-ZST0017
 70-1075-ZST0018	4.5 mmP x 3 mmH Cuff [#8674]	Ø4.7/5.2	70-1075-ZST0018
 70-1075-ZST0019	4.5 mmP x 4 mmH Cuff [#8675]	Ø4.7/5.2	70-1075-ZST0019

LOCATOR is a registered trademark of Zest Anchors, LLC.

LOCATOR® Abutment System

	Description	Code
LOCATOR®		
	Male Processing Package (2-pack) [#8519-2] <i>(Includes: Denture Cap with Black Processing Male, Clear Nylon Male (5lbs), Pink (LR) Less-Retentive Nylon Male (3lbs), Blue Extra-Light-Retention Nylon Male (1.5lb), and White Block-out Spacer.)</i>	70-1075-ZST0024
	Extended Range Male Processing Package (2-pack) [#8540-2] <i>(Includes: Denture Cap with Black Processing Male, Extended-Range Replacement Males: Green (3-4lbs), Orange (2lbs), Red (0.5-1.5lbs), and White Block-Out Spacer.)</i>	70-1075-ZST0025
	Impression Coping (4-pack) [#8505]	70-1075-ZST0007
	Female Analog 4 mmD (4-pack) [#8530]	70-1075-ZST0009
	Core Tool Abutment Holder Sleeve (4-pack) [#8394]	70-1075-ZST0005
	Angle Measurement Guide [#9530]	70-1075-ZST0026
	Processing Spacer (4-pack) [#8569]	70-1075-ZST0010
	Core Tool (Complete) [#8393]	70-1075-ZST0004
	Core Tool Male Removal (Tip End) [#8397]	70-1075-ZST0006
	Parallel Post (4-pack) [#8517]	70-1075-ZST0008
	Core Tool Abutment Driver (Gold End) [#8390]	70-1075-ZST0003
	35 Ncm Torque Wrench Insert Driver (15 mm) [#8260]	70-1075-ZST0001
	35 Ncm Torque Wrench Insert Driver (21 mm) [#8280]	70-1075-ZST0002
	Latch Type Torque Wrench Insert Driver (23 mm) [#8913]	70-1075-ZST0020
	Latch Type Torque Wrench Insert Driver (29 mm) [#8914]	70-1075-ZST0021
	Square Drive Torque Wrench Insert Driver (15 mm) [#8926]	70-1075-ZST0022
	Square Drive Torque Wrench Insert Driver (21 mm) [#8927]	70-1075-ZST0023

LOCATOR is a registered trademark of Zest Anchors, LLC.


Surgical Kit

Description

Code

Inclusive® Tapered Implant Surgical Instrumentation Kit

70-1071-IMP0002



(Bone Profilers sold separately)


Prosthetic Kit

Description

Code

Inclusive® Tapered Implant Prosthetic Kit

70-1072-PRS0001



Features and Benefits

The instrumentation kits allow the clinician to easily organize, store and transport the surgical and prosthetic tooling of the Inclusive Tapered Implant System. The autoclavable kits have been designed for ease of use, with component labels for identification.

The color-coded surgical kit clearly maps the drilling protocol for each implant. Drills are arranged from top to bottom in order of increasing diameter, with color-coded fields to indicate the corresponding diameter of Inclusive Tapered Implant for both soft bone and dense bone protocols. **Please note that Bone Profilers are sold separately (refer to diagram, page 17).**

The prosthetic kit is designed specifically for restorative doctors, organizing all of the necessary prosthetic tooling to help simplify restorative procedures and instrument sterilization.

All instruments associated with the Inclusive Tapered Implant System are machined from corrosion-resistant, surgical stainless steel, and feature standard connectivity.



Sterilization

It is recommended that the instrumentation kits be steam-sterilized for twenty (20) minutes at 132°C/270°F minimum.

General Cleaning Tips:

- Observe universal precautions for the handling of contaminated or biohazardous materials.
- Clean promptly after each use, to prevent biological fluids and tissues from drying on the instruments.
- When applicable, disassemble parts and instruments prior to cleaning.
- Do not rely solely on automatic cleaning. Thorough manual cleaning is recommended.
- Preliminary cleaning should consist of wiping parts, soaking them in a lukewarm enzymatic solution for a minimum of twenty (20) minutes, and rinsing them with running water.
- Routine cleaning should consist of (a) washing parts using a broad spectrum cleaning solution, followed by thorough rinsing and drying; and (b) sonicate parts fully submerged in cleaning solution for at least ten (10) minutes, preferably at 45-50kHz, followed by thorough rinsing and drying.

Drills & Taps

	Description	Code
	Lance Drill Ø1.5 mm	70-1071-SRG0030
	Pilot Drill Ø2.3/2.0 mm, Long	70-1071-SRG0007
	Pilot Drill Ø2.3/2.0 mm, Short	70-1071-SRG0009
	Surgical Drill Ø2.8/2.3 mm, Long	70-1071-SRG0011
	Surgical Drill Ø2.8/2.3 mm, Short	70-1071-SRG0013
	Surgical Drill Ø3.4/2.8 mm, Long	70-1071-SRG0015
	Surgical Drill Ø3.4/2.8 mm, Short	70-1071-SRG0016
	Surgical Drill Ø3.8/3.4 mm, Long	70-1071-SRG0018
	Surgical Drill Ø3.8/3.4 mm, Short	70-1071-SRG0020
	Surgical Drill Ø4.4/3.8 mm, Long	70-1071-SRG0022
	Surgical Drill Ø4.4/3.8 mm, Short	70-1071-SRG0024
	Surgical Drill Ø4.9/4.4 mm, Long	70-1071-SRG0026
	Surgical Drill Ø4.9/4.4 mm, Short	70-1071-SRG0028
	Screw Tap Ø3.2 mm	70-1071-SRG0154
	Screw Tap Ø3.7 mm	70-1071-SRG0031
	Screw Tap Ø4.2 mm	70-1071-SRG0152
	Screw Tap Ø4.7 mm	70-1071-SRG0032
	Screw Tap Ø5.2 mm	70-1071-SRG0033

Drivers & Accessories

	Description	Code
	Inclusive® Bone Profiler Ø3.0 mmP	70-1071-SRG0158
	Bone Profiler Ø3.5 mmP	70-1071-SRG0159
	Bone Profiler Ø4.5 mmP	70-1071-SRG0160
	Tapered Implant Driver Ø3.0 mmP, Long	70-1071-SRG0153
	Tapered Implant Driver Ø3.0 mmP, Short	70-1071-SRG0163
	Tapered Implant Driver Ø3.5/4.5 mmP, Long	70-1071-SRG0161
	Tapered Implant Driver Ø3.5/4.5 mmP, Short	70-1071-SRG0162
	Prosthetic Driver .050" Hex, Long	70-1071-SRG0047
	Prosthetic Driver .050" Hex, Short	70-1071-SRG0048
	Adjustable Torque Wrench	70-1071-SRG0155
	Ratchet Wrench	70-1071-SRG0157
	Parallel Pin	70-1071-SRG0005
	Drill Extender	70-1071-SRG0014
	Driver Handle	70-1071-SRG0049
	Latch-Square Adaptor	70-1072-PRC0050
	Tapered Implant Surgical Adaptor	70-1071-SRG0156
	Tapered Implant Radiographic Template	70-1071-IMP0001


INCLUSIVE®

MINI IMPLANT SYSTEM




Technical Specifications

Below is an enhanced view of the Inclusive Mini Implant with O-ball Head prosthetic attachment, including a side-by-side comparison of the 2.2 mm diameter, 2.5 mm diameter, and 3.0 mm diameter implants. The stated length (10 mm, 13 mm, or 15 mm) refers to the thread length (endosseous section) of the implant, which includes the micro-threaded base of the transgingival collar.



Legacy Mini Implants with O-Ball Head

The original Inclusive Mini Implant design featured a 0.75 mm “short” neck, compared to the 1.5 mm “tall” neck featured in current Inclusive Mini Implant designs. The neck is the instrumentation square located directly above the transgingival collar.



Description	Code
Inclusive® Mini Implant with Collar 2.2 x 10 mm O-Ball	2210
Inclusive® Mini Implant with Collar 2.2 x 13 mm O-Ball	2213
Inclusive® Mini Implant with Collar 2.2 x 15 mm O-Ball	2215
Inclusive® Mini Implant with Collar 2.5 x 10 mm O-Ball	2510
Inclusive® Mini Implant with Collar 2.5 x 13 mm O-Ball	2513
Inclusive® Mini Implant with Collar 2.5 x 15 mm O-Ball	2515
Inclusive® Mini Implant with Collar 3.0 x 10 mm O-Ball	3010
Inclusive® Mini Implant with Collar 3.0 x 13 mm O-Ball	3013
Inclusive® Mini Implant with Collar 3.0 x 15 mm O-Ball	3015


Mini Implants with O-Ball Head

Inclusive Mini Implants are self-tapping threaded titanium screws indicated for long-term applications. Inclusive Mini Implants may also be used for provisional applications. These devices will allow immediate loading and long-term stabilization of dentures and provisional stabilization of dentures while standard implants heal. The implants are to be used for immediate loading only in the presence of primary stability and appropriate occlusal loading.

Description			Code
Inclusive® Mini Implant O-Ball			
 70-1068-IMP0001	 70-1068-IMP0002	 70-1068-IMP0003	2.2 mmD x 10 mmL
			2.2 mmD x 13 mmL
			2.2 mmD x 15 mmL
 70-1068-IMP0004	 70-1068-IMP0005	 70-1068-IMP0006	2.5 mmD x 10 mmL
			2.5 mmD x 13 mmL
			2.5 mmD x 15 mmL
 70-1068-IMP0007	 70-1068-IMP0008	 70-1068-IMP0009	3.0 mmD x 10 mmL
			3.0 mmD x 13 mmL
			3.0 mmD x 15 mmL

Prosthetic Components

The O-ring Housing serves as the retention cap for each Inclusive Mini Implant that is placed. The housing is picked up by the patient's prosthesis using a hard reline of the new or existing denture. These loading procedures can typically be performed immediately after placement of the Inclusive Mini Implants, provided primary stability and appropriate occlusal loading are assured.



The housing is 3.5 mm in height, and 4.75 mm in diameter. The same housing is used regardless of the diameter or length of the selected Inclusive Mini Implant. A vertical space of approximately 1 mm between the crest of the implant collar and the mouth of the housing provides for up to 30 degrees of angular divergence between seated implants.

Each housing comes packaged with one 3/16" x 1/16" O-ring seated within, plus one extra replacement O-ring.

Description	Code
-------------	------

Inclusive® Mini Implant



70-1066-MNP0006

O-Ring Housing	70-1066-MNP0006
----------------	-----------------



70-1066-MNP0008

Prosthetic O-Rings 3/16" x 1/16" (10-pack)	70-1066-MNP0008
--	-----------------



70-1066-MNI0001

Blockout Shims (25-pack)	70-1066-MNI0001
--------------------------	-----------------



70-1066-MNP0004

Impression Coping (4-pack)	70-1066-MNP0004
----------------------------	-----------------



70-1066-MNP0005

Replica (4-pack)	70-1066-MNP0005
------------------	-----------------



70-1066-MNP0007

Instrumentation O-Rings 2.5 mm x 1 mm (5-pack)	70-1066-MNP0007
--	-----------------

Legacy Prosthetic Components

The original Inclusive Mini Implant design featured a 0.75 mm "short" neck, compared to the 1.5 mm "tall" neck featured in current Inclusive Mini Implant designs. The neck is the instrumentation square located directly above the transgingival collar. The components below are for use during the restorative phase of cases utilizing the original "short" neck (0.75 mm) design.

Description	Code
-------------	------

Inclusive® Mini Implant Legacy



70-1066-MNP0003


Impression Coping .75 mm Neck (4-pack)	70-1066-MNP0003
--	-----------------



70-1066-MNP0002

Replica .75 mm Neck (4-pack)	70-1066-MNP0002
------------------------------	-----------------

Surgical Kit Overview



Surgical Kit

All instruments associated with the Inclusive Mini Implant System are machined from high-quality, corrosion-resistant, surgical stainless steel, and feature standard connectivity.

Description	Code
Inclusive® Mini Implant Surgical Instrumentation Kit	70-1067-MNI0006



Features and Benefits

The Mini Implant Surgical Instrumentation Kit allows the clinician to easily organize, store and transport the components of the Inclusive Mini Implant System. It has been designed for ease of use, with component markings and angled rows for easy identification and to help maximize space. The durable, synthetic material is autoclavable and provides a snug, secure fit.

Drills are arranged from left to right in order of increasing diameter, with color-coded fields to indicate the corresponding diameter of Inclusive Mini Implant.

Sterilization

It is recommended that the Surgical Instrumentation Kit be steam-sterilized for twenty (20) minutes at 132°C/270°F minimum.

General Cleaning Tips:

- Observe universal precautions for the handling of contaminated or biohazardous materials.
- Clean promptly after each use, to prevent biological fluids and tissues from drying on the instruments.
- When applicable, disassemble parts and instruments prior to cleaning.
- Do not rely solely on automatic cleaning. Thorough manual cleaning is recommended.
- Preliminary cleaning should consist of wiping parts, soaking them in a lukewarm enzymatic solution for a minimum of twenty (20) minutes, and rinsing them with running water.
- Routine cleaning should consist of (a) washing parts using a broad spectrum cleaning solution, followed by thorough rinsing and drying; and (b) sonicate parts fully submerged in cleaning solution for at least ten (10) minutes, preferably at 45-50 kHz, followed by thorough rinsing and drying.

Instrumentation

	Description	Code
Inclusive® Mini Implant		
	Cortical Bone Drill Ø1.5 mm	70-1071-SRG0008
	Cortical Bone Drill Ø1.7 mm	70-1071-SRG0010
	Cortical Bone Drill Ø2.4 mm	70-1071-SRG0012
	Drill Extender	70-1071-SRG0014
	Mini Implant Driver 13 mm	70-1067-MNI0004
	Mini Implant Driver 5 mm	70-1067-MNI0003
	Handpiece Mini Implant Driver 13 mm	70-1067-MNI0001
	Handpiece Mini Implant Driver 5 mm	70-1067-MNI0002
	Ratchet Wrench	70-1071-SRG0157
	Latch-Square Adaptor	70-1072-PRC0050
	Mini Implant Sterilization Tray & Lid	70-1067-MNI0005
	Mini Implant Radiographic Template	70-1067-MNI0007

Implant Packaging

Inclusive Implants are shipped sterile, and should not be resterilized. They are for single-use only, prior to the expiration date. Do not use implants if packaging has been compromised or previously opened. Do not handle implant surfaces directly. Users are advised to visually inspect packaging to ensure seals and contents are intact prior to use. Please refer to the individual product label for all relevant product information and cautions.

Explanation of Label Codes:

1. Official product description
2. Product image
3. Reference number (product code)
4. Lot number
5. By prescription only
6. Do not use if tampered with
7. Gamma sterilization symbol
8. Expiration date
9. Consult Instructions for Use (IFU)
10. Do not re-sterilize
11. For single-use only
12. Manufacturer
13. Country of origin
14. FDA Unique Device Identification (UDI)
15. Notified body number
16. European Authorized Representative
17. Date of manufacture



1 Inclusive® Tapered Implant
5.2 mmD x 13 mmL x 4.5 mmP

1 pc. **LOT** 1234567 **2**

7 **STERILE R**

6 Contents sterile unless packaging is opened or damaged

12 Prismatic Dentalcraft Inc.
2212 Dupont Drive
Irvine, CA 92612

17 YYYY-MM-DD

8 Use by YYYY-MM-DD

11 Do not re-use

10 Do not re-sterilize

13 Made in USA

5 Rx only

9 Consult instructions for use

15 CE 0086

16 MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany

14

3 **REF** 70-1070-IMP0018

inCLUSIVE®
Tapered Implant
5.2 mmD x 13 mmL x 4.5 mmP

4

15

3011535_6.0

+D745701070IMP00180/SS80132104051234567/16D20160405K

Policies and Warranty

Product Return Policy

Products may be returned at the customer's expense for credit within 30 days of invoice date. All returned products must meet the following conditions:

- A copy of the original invoice must accompany the products.
- Products must be packaged to arrive at the seller's facility undamaged.
- Discontinued, obsolete, expired, damaged, or opened items will not be accepted for return.
- Amount credited will be based on invoice price, less 15 percent for restocking fee.
- Shipping charges are the responsibility of the customer and will not be credited.

Product & Pricing Changes

Because products and equipment are continually undergoing refinement in design and manufacturing methods, we reserve the right to improve, modify, or discontinue products and equipment or change pricing at any time without incurring any obligation and without prior notice.

Warranty

Limited Warranty—Prismatik Dentalcraft, Inc.

Prismatik Dentalcraft, Inc. ("Prismatik"), is the manufacturer of dental products (the "product"), including Inclusive® Dental Implants ("implants"). For a period from the original purchase date of seven (7) years for implants and six (6) months for ceramic blanks and any other product ("the warranty period"), Prismatik will at its option replace or refund the purchase price of any product, to the original purchaser ("user"), that is returned due to defects in material and manufacture.

NO GUARANTEE OR WARRANTY IS IMPLIED OTHER THAN EXPRESSLY STATED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

Prismatik shall not be liable for any incidental or consequential damages, whether foreseeable or not, caused by defects in the product or dental devices produced using said product. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, user's exclusive remedy and Prismatik's sole obligation shall be replacement or refund of the purchase price of the product. For replacement or refund under this warranty, the original purchaser shall send the product at its own expense, postage prepaid, to the seller.

Parts Index

CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE
INCLUSIVE TAPERED IMPLANT SYSTEM					
IMPLANTS					
70-1070-IMP0036.....	Tapered Implant 3.2 mmD x 10 mmL.....	7	70-1071-SRG0156.....	Surgical Adaptor.....	20
70-1070-IMP0032.....	Tapered Implant 3.2 mmD x 11.5 mmL.....	7	70-1071-SRG0153.....	Tapered Implant Driver Ø3.0 mmP, Long.....	20
70-1070-IMP0033.....	Tapered Implant 3.2 mmD x 13 mmL.....	7	70-1071-SRG0163.....	Tapered Implant Driver Ø3.0 mmP, Short.....	20
70-1070-IMP0034.....	Tapered Implant 3.2 mmD x 16 mmL.....	7	70-1071-SRG0161.....	Tapered Implant Driver Ø3.5/4.5 mmP, Long.....	20
70-1070-IMP0035.....	Tapered Implant 3.2 mmD x 8 mmL.....	7	70-1071-SRG0162.....	Tapered Implant Driver Ø3.5/4.5 mmP, Short.....	20
70-1070-IMP0006.....	Tapered Implant 3.7 mmD x 10 mmL.....	7	PROSTHETICS - HEALING/TEMPORIZATION		
70-1070-IMP0007.....	Tapered Implant 3.7 mmD x 11.5 mmL.....	7	70-1069-IMP0009.....	Healing Abutment 3.0 mmP x 3.2 mmD x 3 mmH.....	8
70-1070-IMP0008.....	Tapered Implant 3.7 mmD x 13 mmL.....	7	70-1069-IMP0010.....	Healing Abutment 3.0 mmP x 3.2 mmD x 5 mmH.....	8
70-1070-IMP0009.....	Tapered Implant 3.7 mmD x 16 mmL.....	7	70-1069-IMP0011.....	Healing Abutment 3.0 mmP x 3.2 mmD x 7 mmH.....	8
70-1070-IMP0005.....	Tapered Implant 3.7 mmD x 8 mmL.....	7	70-1069-IMP0012.....	Healing Abutment 3.0 mmP x 3.7 mmD x 3 mmH.....	8
70-1070-IMP0027.....	Tapered Implant 4.2 mmD x 10 mmL.....	7	70-1069-IMP0013.....	Healing Abutment 3.0 mmP x 3.7 mmD x 5 mmH.....	8
70-1070-IMP0028.....	Tapered Implant 4.2 mmD x 11.5 mmL.....	7	70-1069-IMP0014.....	Healing Abutment 3.0 mmP x 3.7 mmD x 7 mmH.....	8
70-1070-IMP0029.....	Tapered Implant 4.2 mmD x 13 mmL.....	7	70-1069-IMP0001.....	Healing Abutment 3.5 mmP x 3.5 mmD x 3 mmH.....	8
70-1070-IMP0030.....	Tapered Implant 4.2 mmD x 16 mmL.....	7	70-1069-IMP0002.....	Healing Abutment 3.5 mmP x 3.5 mmD x 5 mmH.....	8
70-1070-IMP0026.....	Tapered Implant 4.2 mmD x 8 mmL.....	7	70-1069-IMP0015.....	Healing Abutment 3.5 mmP x 3.5 mmD x 7 mmH.....	8
70-1070-IMP0011.....	Tapered Implant 4.7 mmD x 10 mmL.....	7	70-1069-IMP0003.....	Healing Abutment 3.5 mmP x 4.7 mmD x 3 mmH.....	8
70-1070-IMP0012.....	Tapered Implant 4.7 mmD x 11.5 mmL.....	7	70-1069-IMP0004.....	Healing Abutment 3.5 mmP x 4.7 mmD x 5 mmH.....	8
70-1070-IMP0013.....	Tapered Implant 4.7 mmD x 13 mmL.....	7	70-1069-IMP0016.....	Healing Abutment 3.5 mmP x 4.7 mmD x 7 mmH.....	8
70-1070-IMP0014.....	Tapered Implant 4.7 mmD x 16 mmL.....	7	70-1069-IMP0005.....	Healing Abutment 4.5 mmP x 4.5 mmD x 3 mmH.....	8
70-1070-IMP0010.....	Tapered Implant 4.7 mmD x 8 mmL.....	7	70-1069-IMP0006.....	Healing Abutment 4.5 mmP x 4.5 mmD x 5 mmH.....	8
70-1070-IMP0016.....	Tapered Implant 5.2 mmD x 10 mmL.....	7	70-1069-IMP0017.....	Healing Abutment 4.5 mmP x 4.5 mmD x 7 mmH.....	8
70-1070-IMP0017.....	Tapered Implant 5.2 mmD x 11.5 mmL.....	7	70-1069-IMP0007.....	Healing Abutment 4.5 mmP x 5.7 mmD x 3 mmH.....	8
70-1070-IMP0018.....	Tapered Implant 5.2 mmD x 13 mmL.....	7	70-1069-IMP0008.....	Healing Abutment 4.5 mmP x 5.7 mmD x 5 mmH.....	8
70-1070-IMP0019.....	Tapered Implant 5.2 mmD x 16 mmL.....	7	70-1069-IMP0018.....	Healing Abutment 4.5 mmP x 5.7 mmD x 7 mmH.....	8
70-1070-IMP0015.....	Tapered Implant 5.2 mmD x 8 mmL.....	7	70-1073-PRA0048.....	Temporary Abutment 3.0 mmP.....	9
INSTRUMENTATION			70-1073-PRA0022.....	Temporary Abutment 3.5 mmP.....	9
70-1071-PRS0001.....	Prosthetic Kit.....	18	70-1073-PRA0024.....	Temporary Abutment 4.5 mmP.....	9
70-1071-IMP0002.....	Surgical Instrumentation Kit.....	17	70-1073-PRA0047.....	Temporary Abutment, Non-Engaging, 3.0 mmP.....	9
INSTRUMENTATION - DRILLS & TAPS			70-1073-PRA0023.....	Temporary Abutment, Non-Engaging, 3.5 mmP.....	9
70-1071-SRG0030.....	Lance Drill Ø1.5 mm.....	19	70-1073-PRA0025.....	Temporary Abutment, Non-Engaging, 4.5 mmP.....	9
70-1071-SRG0007.....	Pilot Drill Ø2.3/2.0 mm, Long.....	19	PROSTHETICS - IMPRESSIONS/ANALOGS		
70-1071-SRG0009.....	Pilot Drill Ø2.3/2.0 mm, Short.....	19	70-1072-PRC0055.....	Analog 3.0 mmP.....	10
70-1071-SRG0154.....	Screw Tap Ø3.2 mm.....	19	70-1072-PRC0015.....	Analog 3.5 mmP.....	10
70-1071-SRG0031.....	Screw Tap Ø3.7 mm.....	19	70-1072-PRC0017.....	Analog 4.5 mmP.....	10
70-1071-SRG0152.....	Screw Tap Ø4.2 mm.....	19	70-1073-PRA0051.....	Clinical Scanning Abutment, Anterior, 3.0 mmP.....	10
70-1071-SRG0032.....	Screw Tap Ø4.7 mm.....	19	70-1073-PRA0042.....	Clinical Scanning Abutment, Anterior, 3.5 mmP.....	10
70-1071-SRG0033.....	Screw Tap Ø5.2 mm.....	19	70-1073-PRA0043.....	Clinical Scanning Abutment, Anterior, 4.5 mmP.....	10
70-1071-SRG0011.....	Surgical Drill Ø2.8/2.3 mm, Long.....	19	70-1073-PRA0050.....	Clinical Scanning Abutment, Posterior, 3.0 mmP.....	10
70-1071-SRG0013.....	Surgical Drill Ø2.8/2.3 mm, Short.....	19	70-1073-PRA0040.....	Clinical Scanning Abutment, Posterior, 3.5 mmP.....	10
70-1071-SRG0015.....	Surgical Drill Ø3.4/2.8 mm, Long.....	19	70-1073-PRA0041.....	Clinical Scanning Abutment, Posterior, 4.5 mmP.....	10
70-1071-SRG0016.....	Surgical Drill Ø3.4/2.8 mm, Short.....	19	70-1073-PRA0049.....	Laboratory Scanning Abutment 3.0 mmP (3/pk).....	10
70-1071-SRG0018.....	Surgical Drill Ø3.8/3.4 mm, Long.....	19	70-1073-PRA0033.....	Laboratory Scanning Abutment 3.5 mmP (3/pk).....	10
70-1071-SRG0020.....	Surgical Drill Ø3.8/3.4 mm, Short.....	19	70-1073-PRA0034.....	Laboratory Scanning Abutment 4.5 mmP (3/pk).....	10
70-1071-SRG0022.....	Surgical Drill Ø4.4/3.8 mm, Long.....	19	70-1072-PRC0051.....	Transfer Coping, Closed-Tray, 3.0 mmP x 3.2 mmD.....	9
70-1071-SRG0024.....	Surgical Drill Ø4.4/3.8 mm, Short.....	19	70-1072-PRC0053.....	Transfer Coping, Closed-Tray, 3.0 mmP x 3.7 mmD.....	9
70-1071-SRG0026.....	Surgical Drill Ø4.9/4.4 mm, Long.....	19	70-1072-PRC0007.....	Transfer Coping, Closed-Tray, 3.5 mmP x 3.5 mmD.....	9
70-1071-SRG0028.....	Surgical Drill Ø4.9/4.4 mm, Short.....	19	70-1072-PRC0009.....	Transfer Coping, Closed-Tray, 3.5 mmP x 4.7 mmD.....	9
INSTRUMENTATION - DRIVERS & ACCESSORIES			70-1072-PRC0011.....	Transfer Coping, Closed-Tray, 4.5 mmP x 4.5 mmD.....	9
70-1071-SRG0155.....	Adjustable Torque Wrench.....	20	70-1072-PRC0013.....	Transfer Coping, Closed-Tray, 4.5 mmP x 5.7 mmD.....	9
70-1071-SRG0158.....	Bone Profiler Ø3.0.....	20	70-1072-PRC0052.....	Transfer Coping, Open-Tray, 3.0 mmP x 3.2 mmD.....	10
70-1071-SRG0159.....	Bone Profiler Ø3.5.....	20	70-1072-PRC0054.....	Transfer Coping, Open-Tray, 3.0 mmP x 3.7 mmD.....	10
70-1071-SRG0160.....	Bone Profiler Ø4.5.....	20	70-1072-PRC0006.....	Transfer Coping, Open-Tray, 3.5 mmP x 3.5 mmD.....	10
70-1071-SRG0014.....	Drill Extender.....	20,27	70-1072-PRC0008.....	Transfer Coping, Open-Tray, 3.5 mmP x 4.7 mmD.....	10
70-1071-SRG0049.....	Driver Handle.....	20	70-1072-PRC0010.....	Transfer Coping, Open-Tray, 4.5 mmP x 4.5 mmD.....	10
70-1072-PRC0050.....	Latch-Square Adaptor.....	20,27	70-1072-PRC0012.....	Transfer Coping, Open-Tray, 4.5 mmP x 5.7 mmD.....	10
70-1071-SRG0005.....	Parallel Pin.....	20	PROSTHETICS - LOCATOR ABUTMENT SYSTEM		
70-1071-SRG0047.....	Prosthetic Driver .050" Hex, Long.....	20	70-1075-ZST0001.....	35 Ncm Torque Wrench Insert Driver (15 mm).....	16
70-1071-SRG0048.....	Prosthetic Driver .050" Hex, Short.....	20	70-1075-ZST0002.....	35 Ncm Torque Wrench Insert Driver (21 mm).....	16
70-1071-IMP0001.....	Radiographic Template.....	20	70-1075-ZST0026.....	Angle Measurement Guide.....	16
70-1071-SRG0157.....	Ratchet Wrench.....	20,27	70-1075-ZST0004.....	Core Tool (Complete).....	16
			70-1075-ZST0003.....	Core Tool Abutment Driver (Gold End).....	16
			70-1075-ZST0005.....	Core Tool Abutment Holder Sleeve (4/pk).....	16
			70-1075-ZST0006.....	Core Tool Male Removal (Tip End).....	16
			70-1075-ZST0025.....	Extended Range Male Processing Package (2/pk).....	16

Parts Index

CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE
70-1075-ZST0009.....	Female Analog 4 mmD (4/pk)	16	70-1073-PRA0008.....	Titanium Abutment 4.5 mmP x 6 mmH.....	11
70-1075-ZST0007.....	Impression Coping (4/pk)	16	PROSTHETICS - TITANIUM ESTHETIC ABUTMENTS		
70-1075-ZST0020.....	Latch Type Torque Wrench Insert Driver (23 mm)	16	70-1073-PRA0054.....	15° Ti. Esthetic Abut., Ant., 3.0 mmP x 1.5 mmH	11
70-1075-ZST0021.....	Latch Type Torque Wrench Insert Driver (29 mm)	16	70-1073-PRA0059.....	15° Ti. Esthetic Abut., Ant., 3.0 mmP x 2.3 mmH	11
70-1075-ZST0065.....	Locator Abutment 3.0 mmP x 1 mmH Cuff	15	70-1073-PRA0009.....	15° Ti. Esthetic Abut., Anterior, 3.5 mmP	11
70-1075-ZST0066.....	Locator Abutment 3.0 mmP x 3 mmH Cuff	15	70-1073-PRA0031.....	15° Ti. Esthetic Abut., Anterior, 3.5 mmP x 1.5 mmH....	11
70-1075-ZST0011.....	Locator Abutment 3.5 mmP x 0 mmH Cuff	15	70-1073-PRA0010.....	15° Ti. Esthetic Abut., Anterior, 4.5 mmP	11
70-1075-ZST0012.....	Locator Abutment 3.5 mmP x 1 mmH Cuff	15	70-1073-PRA0062.....	15° Ti. Esthetic Abut., Anterior, 4.5 mmP x 1.5 mmH....	11
70-1075-ZST0013.....	Locator Abutment 3.5 mmP x 2.5 mmH Cuff	15	70-1073-PRA0013.....	15° Ti. Esthetic Abut., Posterior, 4.5 mmP.....	11
70-1075-ZST0014.....	Locator Abutment 3.5 mmP x 3.5 mmH Cuff	15	70-1073-PRA0055.....	Ti. Esthetic Abut., Anterior, 3.0 mmP x 1.5 mmH	11
70-1075-ZST0015.....	Locator Abutment 4.5 mmP x 0 mmH Cuff	15	70-1073-PRA0060.....	Ti. Esthetic Abut., Anterior, 3.0 mmP x 2.3 mmH	11
70-1075-ZST0016.....	Locator Abutment 4.5 mmP x 1 mmH Cuff	15	70-1073-PRA0011.....	Ti. Esthetic Abut., Anterior, 3.5 mmP	11
70-1075-ZST0017.....	Locator Abutment 4.5 mmP x 2 mmH Cuff	15	70-1073-PRA0032.....	Ti. Esthetic Abut., Anterior, 3.5 mmP x 1.5 mmH	11
70-1075-ZST0018.....	Locator Abutment 4.5 mmP x 3 mmH Cuff	15	70-1073-PRA0012.....	Ti. Esthetic Abut., Anterior, 4.5 mmP	11
70-1075-ZST0019.....	Locator Abutment 4.5 mmP x 4 mmH Cuff	15	70-1073-PRA0063.....	Ti. Esthetic Abut., Anterior, 4.5 mmP x 1.5 mmH	11
70-1075-ZST0024.....	Male Processing Package (2/pk)	16	70-1073-PRA0014.....	Ti. Esthetic Abut., Posterior, 4.5 mmP	11
70-1075-ZST0008.....	Parallel Post (4/pk)	16	70-1073-PRA0015.....	Ti. Esthetic Abut., Posterior, Wide, 4.5 mmP	11
70-1075-ZST0010.....	Processing Spacer (4/pk)	16	PROSTHETICS - UCLA ABUTMENTS		
70-1075-ZST0022.....	Square Drive Torque Wrench Insert Driver (15 mm)	16	70-1073-PRA0058.....	UCLA Plastic Abutment, 3.0 mmP	12
70-1075-ZST0023.....	Square Drive Torque Wrench Insert Driver (21 mm)	16	70-1073-PRA0018.....	UCLA Plastic Abutment, 3.5 mmP	12
PROSTHETICS - MULTI-UNIT ABUTMENT SYSTEM			70-1073-PRA0020.....	UCLA Plastic Abutment, 4.5 mmP	12
70-1074-PRS0001.....	17° Multi-Unit Abutment 3.5 mmP x 2 mmH	13	70-1073-PRA0057.....	UCLA Plastic Abutment, Non-Engaging, 3.0 mmP	12
70-1074-PRS0002.....	17° Multi-Unit Abutment 3.5 mmP x 3 mmH	13	70-1073-PRA0019.....	UCLA Plastic Abutment, Non-Engaging, 3.5 mmP	12
70-1074-PRS0003.....	17° Multi-Unit Abutment 4.5 mmP x 2 mmH	13	70-1073-PRA0021.....	UCLA Plastic Abutment, Non-Engaging, 4.5 mmP	12
70-1074-PRS0004.....	17° Multi-Unit Abutment 4.5 mmP x 3 mmH	13	INCLUSIVE MINI IMPLANT SYSTEM		
70-1074-PRS0005.....	17° Multi-Unit Abutment 4.5 mmP x 5 mmH	13	IMPLANTS		
70-1074-PRS0006.....	30° Multi-Unit Abutment 3.5 mmP x 3 mmH	13	70-1068-IMP0001.....	Mini Implant O-Ball 2.2 mmD x 10 mmL	23
70-1074-PRS0007.....	30° Multi-Unit Abutment 4.5 mmP x 3 mmH	13	70-1068-IMP0002.....	Mini Implant O-Ball 2.2 mmD x 13 mmL	23
70-1074-PRS0008.....	30° Multi-Unit Abutment 4.5 mmP x 5 mmH	13	70-1068-IMP0003.....	Mini Implant O-Ball 2.2 mmD x 15 mmL	23
70-1074-PRS0031.....	Angled Multi-Unit Abutment Screw	14	70-1068-IMP0004.....	Mini Implant O-Ball 2.5 mmD x 10 mmL	23
70-1074-PRS0040.....	Multi-Unit Abutment 3.0 mmP x 1.5 mmH	13	70-1068-IMP0005.....	Mini Implant O-Ball 2.5 mmD x 13 mmL	23
70-1074-PRS0041.....	Multi-Unit Abutment 3.0 mmP x 2.5 mmH	13	70-1068-IMP0006.....	Mini Implant O-Ball 2.5 mmD x 15 mmL	23
70-1074-PRS0042.....	Multi-Unit Abutment 3.0 mmP x 3.5 mmH	13	70-1068-IMP0007.....	Mini Implant O-Ball 3.0 mmD x 10 mmL	23
70-1074-PRS0009.....	Multi-Unit Abutment 3.5 mmP x 1 mmH	13	70-1068-IMP0008.....	Mini Implant O-Ball 3.0 mmD x 13 mmL	23
70-1074-PRS0010.....	Multi-Unit Abutment 3.5 mmP x 2 mmH	13	70-1068-IMP0009.....	Mini Implant O-Ball 3.0 mmD x 15 mmL	23
70-1074-PRS0011.....	Multi-Unit Abutment 4.5 mmP x 1 mmH	13	INSTRUMENTATION		
70-1074-PRS0012.....	Multi-Unit Abutment 4.5 mmP x 2 mmH	13	70-1071-SRG0008.....	Cortical Bone Drill Ø1.5 mm.....	27
70-1074-PRS0013.....	Multi-Unit Abutment 4.5 mmP x 3 mmH	13	70-1071-SRG0010.....	Cortical Bone Drill Ø1.7 mm	27
70-1074-PRS0014.....	Multi-Unit Abutment 4.5 mmP x 5 mmH	13	70-1071-SRG0012.....	Cortical Bone Drill Ø2.4 mm.....	27
70-1074-PRS0016.....	Multi-Unit Abutment Analog.....	14	70-1067-MNI0001.....	Handpiece Mini Implant Driver 13 mm.....	27
70-1073-PRA0037.....	Multi-Unit Coping 4.5 mmH, Non-Engaging	14	70-1067-MNI0002.....	Handpiece Mini Implant Driver 5 mm	27
70-1074-PRS0028.....	Multi-Unit Driver.....	14	70-1067-MNI0004.....	Mini Implant Driver 13 mm.....	27
70-1074-PRS0018.....	Multi-Unit Gold Bar Coping.....	14	70-1067-MNI0003.....	Mini Implant Driver 5 mm.....	27
70-1074-PRS0030.....	Multi-Unit Gold/Plastic Coping.....	14	70-1067-MNI0007.....	Radiographic Template	27
70-1074-PRS0061.....	Multi-Unit Guide Pin .050" Hex	14	70-1067-MNI0005.....	Sterilization Tray & Lid	27
70-1074-PRS0021.....	Multi-Unit Impression Coping, Closed-Tray.....	14	70-1067-MNI0006.....	Surgical Instrumentation Kit.....	26
70-1074-PRS0052.....	Multi-Unit Impression Coping, Open-Tray.....	14	PROSTHETICS		
70-1072-PRC0014.....	Multi-Unit Prosthetic Screw .050" Hex.....	14	70-1066-MNI0001.....	Blockout Shims (25/pk).....	24
70-1074-PRS0015.....	Multi-Unit Temporary Healing Cap	14	70-1066-MNP0004.....	Impression Coping (4/pk)	24
70-1074-PRS0025.....	Multi-Unit Titanium Temporary	14	70-1066-MNP0007.....	Instrumentation O-Rings 2.5 mm x 1 mm (5/pk)	24
PROSTHETICS - SCREWS & GUIDE PINS			70-1066-MNP0005.....	Mini Implant Replica (4/pk).....	24
70-1072-PRC0002.....	Guide Pin (5/pk).....	12	70-1066-MNP0006.....	O-Ring Housing	24
70-1072-PRC0056.....	Guide Pin 3.0 mmP (5/pk).....	12	70-1066-MNP0008.....	Prosthetic O-Rings 3/16" x 1/16" (10/pk).....	24
70-1072-PRC0001.....	Titanium Screw	12	PROSTHETICS - LEGACY		
70-1072-PRC0003.....	Titanium Screw (5/pk)	12	70-1066-MNP0003.....	Legacy Impression Coping .75 mm Neck (4/pk)	24
70-1072-PRC0058.....	Titanium Screw 3.0 mmP	12	70-1066-MNP0002.....	Legacy Mini Implant Replica .75 mm Neck (4/pk)	24
70-1072-PRC0057.....	Titanium Screw 3.0 mmP (5/pk)	12			
PROSTHETICS - TITANIUM ABUTMENTS					
70-1073-PRA0052.....	Titanium Abutment 3.0 mmP x 4.5 mmH.....	11			
70-1073-PRA0053.....	Titanium Abutment 3.0 mmP x 6 mmH.....	11			
70-1073-PRA0045.....	Titanium Abutment 3.5 mmP x 4.5 mmH.....	11			
70-1073-PRA0006.....	Titanium Abutment 3.5 mmP x 6 mmH.....	11			
70-1073-PRA0046.....	Titanium Abutment 4.5 mmP x 4.5 mmH.....	11			

INCLUSIVE®

DENTAL IMPLANT SYSTEM



Designed & Manufactured in the USA
by



PRISMATIK
DENTALCRAFT, INC.

(A wholly owned subsidiary of Glidewell Laboratories)

2212 Dupont Dr. • Irvine, CA 92612

CE 0086

EC	REP
----	-----

MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

MKT 787_5.0 (3019959)
REV ID 7175

GD-809-112216