

# SIX LEVELS OF PROTECTION

## JUNIOR



One layer of EVA material (3 mm) with added incisal and occlusal protection. Designed specifically for children with mixed dentition.

## LIGHT



Two layers of laminated EVA (3 mm) material. Designed for wrestling, volleyball, mountain biking and motocross. (Only available in clear, bright yellow, bright red, light blue and green.)

## LIGHT PRO




Three layers of laminated EVA (3.5 mm) material (one hardened layer). Specially designed to provide maximum comfort and protection for all sports.

## MEDIUM



Two layers of laminated EVA (5 mm) material. Specially fabricated for soccer, rugby, basketball, softball, rollerblading and skateboarding.

## HEAVY



Two layers of laminated EVA (5 mm) material with three unique power dispersion bands. Specially designed for baseball, football, racquetball, martial arts and boxing.

## HEAVY PRO



Three layers of laminated EVA (5 mm) material (one hardened layer). Custom fabricated for ice, field and street hockey, kickboxing, and other heavy contact sports where blows from pointed objects are expected.

Helmet straps available in red and black

Personalize **PlaySafe**® mouthguards with one of these standard colors ...

PURE WHITE	BRIGHT YELLOW	GOLD	BRIGHT GREEN
DEEP GREEN	SILVER	LIGHT BLUE	BRIGHT BLUE
DARK BLUE	BRIGHT PINK	BRIGHT RED	DEEP RED
MAROON	DEEP BLACK	◆ RED, WHITE & BLUE	◆ LILAC
◆ TANGERINE	◆ YOLK YELLOW	CAMOUFLAGE	CAMOUFLAGE STRIP
LAVA	LAVA STRIP	TIE-DYE	RAINBOW
CONFETTI	ZEBRA	GOLD FLAKES	SILVER FLAKES
SPORTS	ANTS	TEAM NAMES AND MANY OTHER STICKERS AVAILABLE	
		◆ NOT AVAILABLE FOR LIGHT OR LIGHT PRO	

... or create a custom design!



PlaySafe is a registered trademark of Glidewell Laboratories.



MKT 9977\_1.0  
REV ID 11569  
GL-1129-042018

# PREVENT INJURY

Custom-made, dentist-prescribed sports mouthguard



**PLAYSAFE**®  
sports mouthguard

## WHY CUSTOM FABRICATED SPORTS MOUTHGUARDS?

Laboratory-made PlaySafe® sports mouthguards are prescribed for both professional and amateur athletes because they offer the highest level of protection during sports activity. When worn properly, they can prevent dental trauma, such as fractured teeth. They are also indicated for patients with an orthodontic situation, such as braces, who must wear a mouthguard for preventive measures.

## A CUSTOM FIT ENSURES PROTECTION

PlaySafe sports mouthguards provide superior fit and retention compared to boil-and-bite mouthguards. Available in six levels of protection, the mouthguards range from one to three layers of laminated ethylvinylacetate material with a final occlusal thickness between 3 mm and 5 mm. PlaySafe mouthguards can be customized with stickers and team logos, and are available in a variety of colors to fit the personality of every star athlete in your practice!



PLAYSAFE SPORTS MOUTHGUARD

VS.

BOIL-AND-BITE MOUTHGUARD



### This could have been prevented!

This serious trauma was caused by an elbow blow. It could have been avoided if the basketball player had been wearing a PlaySafe mouthguard.



### PLAYSAFE® TRI-LAMINATION

Three times the protection

