



New West Dental Lab

2033 Swanson Ave. ■ Lake Havasu City, AZ 86403
800-321-1614 ■ Fax 855-250-2830
newwestlab.com

SPECIAL BRUXZIR ZIRCONIA RX OFFER*

- Carefully package your case, including this Rx, and tape box securely closed.
- To schedule shipping pickup, call us at **800-321-1614**.
- Send this Rx with your next BruxZir case to save \$30 per BruxZir unit (limit 3 units).

Offer expires Dec. 31, 2025

**Limit one specially priced Rx per case. Web offer can only be used a maximum of four times per account. Special discount not valid with any other offer.*

Dr. Name _____ Acct. # _____

Patient ID/Name _____ Age _____
First Last

Deliver by 5 p.m. on _____

Enclosed with case: Impressions Models Bite Photos Other: _____

WEB Rx

- BruxZir Full-Strength** (> 1,000 MPa)
- BruxZir Esthetic (870 MPa)

NOTE: Please send a study model on all work involving anterior teeth.

\$30 off BruxZir Full-Strength



* U T D P 6 5 *

\$30 off BruxZir Esthetic



* D J L A 4 8 *

SHADE INFORMATION

Tooth No. _____


Stamp Shade _____

Final Shade _____

OCCUSAL STAINING

None Light** Medium Dark

PONTIC DESIGN



Signature _____ License # _____ Date _____

Submission of this Rx constitutes agreement with limited warranty terms and conditions. See reverse for details.

TERMS AND WARRANTY INFORMATION



All Restorations Made in the USA

We honor VISA, MASTERCARD, AMEX and DISCOVER.

TERMS: Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. **Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.** Prices subject to change without notice. Rx must be enclosed with original case submission.

NO-FAULT REMAKE POLICY: New West is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.


LIMITED WARRANTY/LIMITATION OF LIABILITY. For warranty terms and conditions and limitation of liability, visit newwestlab.com/warranty.



• BruxZir® Restorations

All rush cases must be prescheduled by calling 800-321-1614 before the case is shipped. Time of pickup and delivery may affect turnaround time.

PREPARATION GUIDELINES



BruxZir Esthetic

- A. 1.25 mm ideal reduction (0.7 mm minimum)
- B. Chamfer or modified shoulder margins preferred
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

BruxZir Full-Strength

- A. 1.0 mm ideal reduction (0.5 mm minimum)
- B. Chamfer or shoulder margins preferred. Feather-edge OK
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins